




Tools Needed

T15 torx

5mm Allen key to
undo the caliper

2x Bleed nipple
tools (Tekro only
tool)

Syringe




Undo the bleed nipple bolt with
the T15



Insert bleed nipple into the lever reservoir.



Attach catch pot.



Fill a third of the syringe with mineral oil.




Remove T15 bleed port bolt.



Insert second Tektro bleed nipple.




Attach syringe of oil to the calliper bleed port.



the catch pot.

Push excess air through the system into



Remove caliper bleed nipple & insert screw.



Remove catchpot and bleed nipple & pump lever until free of small bubbles. (top up accordingly)



When you feel the lever has firmed up, insert the bleed port screw and wipe down brake with isopropyl alcohol.